



NOTICE OF AWARD

Date: December 10, 2024

MS. ANNALYN SITYAR-BERNABE

Proprietress/Admin and Finance Manager

CAS Color Application Specialist

Unit 1213 Entrata Tower 1 Urban Complex Condominium 2609, Civic Drive Filinvest City, Alabang, Muntinlupa City

Dear Sir/Madam,

We are pleased to notify you that PhilRice awarded the procurement of the following:

PARTICULARS	DEPT/OFFICE	TOTAL AMOUNT
Supply and Delivery of Laboratory Chemicals, Supplies, and Materials for Corporate Project (4th Quarter) under PB 24-03-33	Plant Breeding and Biotechnology Division (PBBD)	PhP589,780.00

In this connection, you are hereby required to immediately post a performance bond, the amount of which is in accordance with the provision of RA 9184 which was approved by the President of the Philippines on January 10, 2003 as indicated below:

Form of Performance Security	Amount of Performance Security (Equal to the Percentage of the Total Contract Price)
a) Cash or cashier's/manager's check issued by a Universal or Commercial Bank	Five percent (5%)
b) Bank draft/guarantee or irrevocable letter of credit issued by a Universal or Commercial Bank: Provided, however that it shall be confirmed or authenticated by a Universal or Commercial Bank, if issued by a foreign bank.	
c) Surety bond callable upon demand issued by a surety or insurance company duly certified by the Insurance Commission as authorized to issue such security and/or	Thirty percent (30%)
d) Any combination of the foregoing.	Proportionate to share of form with respect to total amount of security

Failure to post the performance bond within ten (10) calendar days upon receipt of this notice would mean a refusal to accept the award. Such is sufficient ground for PhilRice to impose administrative sanctions on your company. The contract/Purchase Order (PO) shall be released upon compliance with the requirement set forth herein. Kindly coordinate with the BAC Secretariat, at telephone number (044) 456-0650 regarding the above mentioned requirement or any inquiries.

Thank you very much.



Better Rice Communities.

Email: prri.mail@mail.philrice.gov.ph Text Center: (+63) 917-111-7423

Websites: www.philrice.gov.ph, www.pinoyrice.com

Social Media: [web.facebook.com/DAPhilRice](https://www.facebook.com/DAPhilRice)

Regional Office: 3rd Floor, ATI building, Elliptical Road, Diliman, Quezon City

Mobile No.: (+63)928-915-9628

NOA Control No.: 2024-11-0767 (C)



Philippine Rice Research Institute

Central Experiment Station

Maligaya, Science City of Muñoz, 3119 Nueva Ecija

Very truly yours,

EDUARDO JIMMY P. GUILANG
Deputy Executive Director for
Research

Conforme:

FERNAN G. ZAPE

SIGNATURE OVER PRINTED NAME

CAS Color Application Specialist

Letter Rice Communities.

Email: prri.mail@mail.philrice.gov.ph Text Center: (+63) 917-111-7423

Websites: www.philrice.gov.ph, www.pinoyrice.com

Social Media: [web.facebook.com/DAPhilRice](https://www.facebook.com/DAPhilRice)

Regional Office: 3rd Floor, ATI building, Elliptical Road, Diliman, Quezon City

Mobile No.: (+63)928-915-9628

NOA Control No.: 2024-11-0767 (C)





List of Goods/Services for Award to CAS Color Application Specialist

ITEM NO.	QTY	UNIT	ITEM DESCRIPTION	UNIT PRICE	TOTAL PRICE
2	100	piece	Semi-Expendable Technical and Scientific Equipment, Crucible Cole-Parmer Wide-Form Crucible with Cover, Porcelain, 20 mL; Cole-Parmer – Item # EW-63100-03 Withstand temperatures up to 1150 °C (2102 °F) Resistance to acids and alkalies Autoclavable 1 year warranty Delivery Period: 30 calendar days	698.00	69,800.00
5	1	pack	Membrane, Nitrocellulose Blotting Membrane Cytiva Cat. no. 10600124 Amersham Protran 0.45 NC nitrocellulose Western blotting membranes Amersham Protran 0.45?m NC 102mm×133mm 10 sheets/PK Delivery Period: 45 Calendar days	14,800.00	14,800.00
21	1	bottle	Molecular biology reagent, Fetal bovine serum HyClone Cat. no. SV30160.03HI 29181940 Research Grade Fetal Bovine Serum (FBS) South American Origin Heat Inactivated 500ML Delivery Period: 60 calendar days	34,800.00	34,800.00
32	4	unit	Semi-Expendable Technical and Scientific Equipment, Pipette Eppendorf Research® plus, 20ul to 200uL Cat. No. 3123000055 Research plus is highly ergonomic and it features very low weight and operation forces, is extremely reliable, secure volume adjustment and fully autoclavable and easy to use. It has also spring-loaded tip cone, a temporary adjustment option, and improved volume display. ACT label: yes Open system (third party consumables) yes Product type: Pipettes Pipetting type: Air-cushion system Number of channels: 1-channel Volume range: 20 – 200 ?L Color code: yellow Volume increment : 0.2 ?L Operating mode: mechanical Delivery content Pipette, epT.I.P.S.® Box 2.0 pipette tips, Operating manual, Certificate, Adjustment tool (Allen key with blue handle), Adjustment seal for temporary adjustment (x5), Grease for pipettes, Locking ring (to prevent spring action)	24,670.00	98,680.00

Better Rice Communities.


Email: prri.mail@mail.philrice.gov.ph Text Center: (+63) 917-111-7423

Websites: www.philrice.gov.ph, www.pinoyrice.com

Social Media: web.facebook.com/DAPhilRice

Liaison Office: 3rd Floor, ATI building, Elliptical Road, Diliman, Quezon City

Mobile No.: (+63)928-915-9628





Philippine Rice Research Institute

Central Experiment Station

Maligaya, Science City of Muñoz, 3119 Nueva Ecija

ITEM NO.	QTY	UNIT	ITEM DESCRIPTION	UNIT PRICE	TOTAL PRICE
			Delivery Period: 60 calendar days		
33	4	unit	Semi-Expendable Technical and Scientific Equipment, Pipette Eppendorf Research® plus Catalog No. 3123000020 1-channel 0.5 – 10 ?L, medium gray, variable volume incl. epT.I.P.S.® Box 2.0 with 96 pipette tips Research plus is highly ergonomic and has secure volume adjustment. It features very low weight and operation forces, is extremely reliable, fully autoclavable and easy to use. It has also spring-loaded tip cone, a temporary adjustment option, and improved volume display. Pipetting type Air-cushion system Number of channels 1-channel Volume range 0.5 – 10 ?L Color code medium gray Volume increment 0.01 ?L Delivery content Pipette, epT.I.P.S.® Box 2.0 pipette tips, Operating manual, Certificate, Adjustment tool (Allen key with blue handle), Adjustment seal for temporary adjustment (x5), Grease for pipettes, Locking ring (to prevent spring action) Delivery Period: 60 calendar days	24,670.00	98,680.00
34	2	unit	Semi-Expendable Technical and Scientific Equipment, Pipette Cole-Parmer 8 Channel Mechanical Adjustable-Volume Pipette, 0.5 to 10 ?l Cat. No. EW-24501-19 eatures click-stop volume settings for simple and quick use Ergonomic, lightweight push button design makes them easy to operate forward and reverse pipetting modes Drop of the pipette is color-coded by volume for easy identification Volume displayed through handgrip window for quick glance convenience Volume: 0.5 to 10 ?L 3 year warranty Specifications & Description AutoclavableYes; Max Volume (?L)10; Min Volume (?L)0.5; Increments (?L)0.1Accuracy (%)±1.5Number Of Channels8 Tip EjectorYesColor Channel Mechanical Adjustable-Volume Pipette	34,490.00	68,980.00

Better Rice Communities.


Email: prri.mail@mail.philrice.gov.ph Text Center: (+63) 917-111-7423

Websites: www.philrice.gov.ph, www.pinoyrice.com

Social Media: [web.facebook.com/DAPHilRice](https://www.facebook.com/DAPHilRice)

Liaison Office: 3rd Floor, ATI building, Elliptical Road, Diliman, Quezon City

Mobile No.: (+63)928-915-9628





Philippine Rice Research Institute

Central Experiment Station

Maligaya, Science City of Muñoz, 3119 Nueva Ecija

ITEM NO.	QTY	UNIT	ITEM DESCRIPTION	UNIT PRICE	TOTAL PRICE
			Delivery content Pipette, epT.I.P.S.® Box 2.0 pipette tips , Operating manual, Certificate, Adjustment tool (Allen key with blue handle), Adjustment seal for temporary adjustment (x5), Grease for pipettes, Locking ring (to prevent spring action) Delivery Period: 60 calendar days		
35	4	unit	Semi-Expendable Technical and Scientific Equipment, Pipette Eppendorf Research® plus2-20uL Cat. no. - 3123000039 Research plus is highly ergonomic and has secure volume adjustment. It features very low weight and operation forces, is extremely reliable, fully autoclavable and easy to use. It has also spring-loaded tip cone, a temporary adjustment option, and improved volume display ACT label: yes Product type: Pipettes Pipetting type: Air-cushion system Number of channels: 1-channel Volume range: 2 – 20 ?L Color code: yellow Volume increment: 0.02 ?L Operating mode: mechanical Spring-loaded tip cone: yes Open system (third party consumables) yes Delivery content Pipette, epT.I.P.S.® Box 2.0 pipette tips , Operating manual, Certificate, Adjustment tool (Allen key with blue handle), Adjustment seal for temporary adjustment (x5), Grease for pipettes, Locking ring (to prevent spring action) Eppendorf Research® plus, 100-1000uL Delivery Period: 60 calendar days	24,670.00	98,680.00
36	4	unit	Semi-Expendable Technical and Scientific Equipment, Pipette Eppendorf Research® plus, 100-1000uL ACT label of sustainability certification for Eppendorf Research® plus pipette, 100-1000 ?L (blue), single-channel, variable Research plus is highly ergonomic and has secure volume adjustment. It features very low weight and operation forces, is extremely reliable, fully autoclavable and easy to use. It has also spring-loaded tip cone, a temporary adjustment option, and improved volume display. Specifications Eppendorf Research ® plus - 3123000063 ACT label yes	24,670.00	98,680.00

Better Rice Communities.


Email: prri.mail@mail.philrice.gov.ph Text Center: (+63) 917-111-7423

Websites: www.philrice.gov.ph, www.pinoyrice.com

Social Media: [web.facebook.com/DAPHilRice](https://www.facebook.com/DAPHilRice)

Liaison Office: 3rd Floor, ATI building, Elliptical Road, Diliman, Quezon City

Mobile No.: (+63)928-915-9628





Philippine Rice Research Institute

Central Experiment Station

Maligaya, Science City of Muñoz, 3119 Nueva Ecija

ITEM NO.	QTY	UNIT	ITEM DESCRIPTION	UNIT PRICE	TOTAL PRICE
			Product type Pipettes Pipetting type Air-cushion system Number of channels 1-channel Volume range 100 – 1,000 ?L Color code blue Volume increment 1 ?L Operating mode mechanical Spring-loaded tip cone yes Open system (third party consumables) yes Delivery content Pipette, epT.I.P.S.® Box 2.0 pipette tips Operating manual, Certificate, Adjustment tool (Allen key with blue handle), Adjustment seal for temporary adjustment Delivery Period: 60 calendar days		
40	20	pack	Tips, Pipette tips GREINER BIO-ONE Sapphire Pipette Tip 1000 ?l, natural, non-sterile, 1000 pcs./bag Item No.: 777350 Description: Sapphire Pipette Tip, 10 ?l, natural, non-sterile, 1000 pcs./bag Sterile: non-sterile Tips are fully autoclavable; non-pyrogenic certified and RNase and DNase free Delivery Period: 30 calendar days	334.00	6,680.00
					589,780.00

Better Rice Communities.

Email: prri.mail@mail.philrice.gov.ph Text Center: (+63) 917-111-7423

Websites: www.philrice.gov.ph, www.pinoyrice.com

Social Media: [web.facebook.com/DAPhilRice](https://www.facebook.com/DAPhilRice)

Liaison Office: 3rd Floor, ATI building, Elliptical Road, Diliman, Quezon City

Mobile No.: (+63)928-915-9628

